

Accessories for individual requirements



The separately removable armrests with ball and socket joint can be adjusted in the front part and support the patient's arm in a variety of positions.



For your patient's safety four guard rails can be attached quickly and easily to the DreamLine.



The held-up support assists your patient when sitting up and provides therefore additional flexibility.



Also the infusion pole can be mounted space saving directly to the back segment of the DreamLine.



The LED lamp offers pleasant light for your patient for reading or working during the treatment.

Standard colour Chrom grey

Amazonas

Feuer

Lima

Atoll

Other colours are available upon request. Please contact us!

Technical Data:

Dimensions

Chair width, without armrests	90 cm
Total length	230 cm
Mattress width	75 cm
Mattress length	215 cm
Total weight	123 kg

Adjustment features

Back segment angle	0° – 75°
Seat segment angle	0° – 25°
Leg segment angle	0° – 25°
Footrest	23 cm
Height adjustment, incl. mattress	55 – 95 cm

Electrical specifications

Protection class	2
Protection category	IPX 4
In compliance with	EN 60601
Motors	5

Safety

Safe maximum load	200 kg
Angle of shock position	- 12°
Time for shock position	20 sec

The Therapy Chair with one piece lying surface – experience more flexibility and comfort



Bionic Medizintechnik GmbH ... a JMS-corporation · Max-Planck-Straße 21 · D-61381 Friedrichsdorf/Germany
Phone: +49(0)6172-75760 · Fax: +49(0)6172-757610 · info@bionic-jms.com · www.bionic-jms.com

Right for technical modifications reserved
B-MOBDLE/02.2011-09



The DreamLine – offers your patients a variety of sitting positions as well as a very relaxing reclining position



The very low entry height facilitates easy entry – also from a wheelchair.



During the treatment the patient can adjust the DreamLine to achieve a pleasant and relaxing comfort position.



For optimum relaxation the DreamLine offers a completely horizontal position.



The shock position is activated by foot pedal, allowing the medical team hands free operation.

The DreamLine further enhances the advantages of a therapy chair by providing a one piece lying surface with more flexibility in its functions. Besides your patient's requirements for a comfortable sitting position - it also provides a very relaxing lying position.

Five motors ensure a smooth and continuous adjustment by an easy to operate remote control, to every desired sitting or lying position. The height can be adjusted electrically between 55 and 95 cm. For an even more comfortable position, depending on the patient's height, there is an electrically adjustable footrest.

Due to the stable base frame the DreamLine is very safe in all possible positions. The equipotential bonding according EN 60601 provides additional safety for your patient.

Underneath the easy to remove mattress there are plastic plates which are resistant to

disinfectants and easy to clean. They provide together with the mattress the unique lying and sitting comfort of the DreamLine.

An adequate area at the base frame offers enough space for accessories or your patient's belongings.

In case of an emergency the shock position can be activated fast and easily via foot pedal. Additionally, the headboard can be removed quickly to have free access to the patient.

Individual accessories such as armrests, infusion pole and held-up support increase the functionality and flexibility of the DreamLine and can be mounted space saving. If your patient would like to read or work during treatment the LED lamp provides a pleasant light. For more safety four easy to attach guard rails are available.



The easy to operate remote control with pre-programmed sitting and lying position ensures a simple adjustment of the individual positions.



With the help of the electrically adjustable footrest the DreamLine can be adapted easily and flexibly to your patient's height.



The two height adjustable columns guarantee high stability in all positions. An adequate size of storage space is provided at the base frame.



The foot pedal to activate the shock position and the central brake are located at the base frame.



The headboard can be removed fast and easily in case of an emergency to allow full access to the patient's head.



Underneath the mattress easy to clean plastic plates provide a stable, safe and hygienic surface.